

When Seconds Count.

WATER-JEL
Fire Blankets are
packaged in a more
contemporary and
ergonomic version of
the "orange canister"
that has become an
icon in Emergency
Burn Treatment.





Modern technology has advanced the engineering and manufacturing process so that a more efficient unit can be manufactured. This allows for quicker response time in an emergency situation.

Available in 3 sizes:

- · 36" x 30"
- 6' x 5'
- · 8' x 6'

In a burn emergency, every second is critical. By saving precious time, you could be saving a life. Water-Jel Fire Blankets are 100% worsted wool, with an intercellular weave capable of absorbing up to 13 times it own weight of the scientifically formulated cooling gel in which it is soaked. This biodegradable, bacteriostatic and water-soluble gel immediately upon contact, begins to extinguish the flame, cool the burn and help stop burn progression and pain. It stabilizes the patient and helps prepare for transport to a burn center or ER.

- These blankets are vital in areas where high heat flames, arc flash, hot liquid, steam or chemicals exist.
- Blankets can also be used to put out fires and protect a rescuer from heat and flames.
- Medical Devices contain no drugs or other active ingredients.
- Used by EMS, Firefighters, Utility companies, Industrial Institutions and the Military around the world.

Improved Canister Design:

- Swift removal of blanket means faster response to a burn emergency
- Black clips have been replaced no more bending, cracking or resistance of old style.
- New Grip & Turn allows simple, quick opening.
- · Easy lock mechanism keeps blanket securely inside
- · Easy grasp handle on larger units makes carrying convenient
- · All physical specs remain unchanged (weight, size, construction integrity)









WATER-JEL FIRE BLANKETS IN CANISTERS AND POUCHES —

Mounting Brackets also available to hang blanket canisters on wall



WATER-JEL HEAT SHIELD Canister

Blanket Dimensions: 8' x 6'

9672-04 Case Pack: 4 Canisters individually boxed

Case Size: 18.5" x 15.75" x 23.75"

Case Weight: 113 lbs. Case Cube: 4.0 cu. ft. Item UPC: 7 43537 96721 6



HEAVY DUTY MOUNTING BRACKET Holds 9672 Heat Shield Canister

TL-10 **Bracket individually packaged**

Weight: 5.8 lbs.



WATER-JEL FIRE BLANKET-PLUS Canister

Blanket Dimensions: 6' x 5'

7260-04 Case Pack: 4 Canisters individually boxed

Case Size: 18.5" x 15.75" x 18"

Case Weight: 82 lbs. Case Cube: 3.04 cu. ft. Item UPC: 7 43537 72601 1



HEAVY DUTY MOUNTING BRACKET Holds 7260 Fire Blanket Plus Canister

TM-10 **Bracket individually packaged**

Weight: 5.8 lbs.



WATER-JEL BURN WRAP Canister

Blanket Dimensions: 3' x 2.5'

3630-04 Case Pack: 4 Canisters individually boxed

Case Size: 15.5" x 15.5" x 6.75"

Case Weight: 26.5 lbs. Case Cube: .94 cu. ft Item UPC: 7 43537 36301 8



HEAVY DUTY MOUNTING BRACKET

Holds 3630 Burn Wrap Canister

TS-10 Bracket individually packaged Weight: 3.4 lbs.



WATER-JEL FIRE BLANKET-PLUS **Pouch**

P7260-04 Case Pack: 4 Pouches individually boxed. Blanket Dimensions: 6' x 5'

Case Size: 19.5" x 11.75" x 11" Case Weight: 56 lbs.

Case Cube: 1.46 cu. ft. Item UPC: 7 43537 07260 6



WATER-JEL BURN WRAP **Pouch**

P3630-4 Case Pack: 4 Pouches

individually boxed. Blanket Dimensions: 3' x 2.5' Case Size: 16" X 9" X 7.75"

Case Weight: 19 lbs. Case Cube: .65 cu. ft.

Item UPC: 7 43537 03630 1

Shelf life on Fire Blankets: 5 years from date of manufacture

Blankets also available in canisters designed to Military/Government spec see your Sales Rep for information



50 Broad Street, Carlstadt, NJ 07072 Tel: (201) 507-8300 • Fax (201) 507-8325 • (800) 275-3433 E-mail: info@waterjel.com • Website: www.waterjel.com